

DGspin

Digital centrifuge for processing DG Gel Cards

Manual Blood Bank System



Proven performance

- Gel agglutination technique provides proven accuracy and reliability
- Manual system utilizes same gel cards and reagents as Grifols fully automated systems

Intelligent design

- Capacity to centrifuge up to 24 gel cards
- Display indicates centrifuge speed and time remaining
- Pre-set centrifuge time of 9 minutes


TYPING

GRIFOLS

KEMENKES RI AKL XXXXXXXX

Distributed by:
PT. Frismed Hoslab Indonesia
 Ruko Bidex Blok G No. 30-31 CBD
 BSD City, Tangerang Selatan - Indonesia 15321
 Phone/Fax : +6221 53167150-52/+6221 5316715

Manufacturer :
Diagnostic Grifols, S.A.
 Pg. Fluvial, 24
 08015 Parets del Vallès
 Barcelona, Spain
 Tel: +34 935 710 400



DGspin

Manual Blood Bank System

Digital centrifuge for processing DG Gel Cards



Technical specifications

SPECIFICATIONS

Height	180 mm
Width	410 mm
Depth	525 mm
Weight	Approx. 9 kg

SAMPLE STATION

Load capacity	24 cards
DG Spin head	Removable and interchangeable
Centrifugation conditions	Pre-programmed for DG Gel cards
Centrifugation speed	990 rpm \pm 10 rpm
Centrifugation speed units displayed	Programmable, rpm or g
Centrifugation time	9 minutes
Alarm	Yes, to indicate end of centrifugation and error situations
Volume of alarm	Programmable, 3 levels
Safety systems	Out-of-balance control / Detection of door opening / Safety closure
Interface	LCD screen and buttons

ELECTRICAL REQUIREMENTS

Voltage	100-240V
Frequency	50-60 Hz
Maximum input power consumption	50W
Fuses	2 x T 2A L 250 V, 5 x 20 mm

ENVIRONMENTAL CONDITIONS

Temperature	15°C - 40°C (59°F - 104°F)
Maximum relative humidity non-condensing	80%
Maximum altitude	3000m
Power supply maximum voltage fluctuations	\pm 10% of the nominal voltage


For more information, ask your local Grifols representative.



KEMENKES RI AKL XXXXXXXX



Distributed by:
PT. Frismed Hoslab Indonesia
Ruko Bidex Blok G No. 30-31 CBD
BSD City, Tangerang Selatan - Indonesia 15321
Phone/Fax : +6221 53167150-52/+6221 5316715

Manufacturer :
Diagnostic Grifols, S.A.
Pg. Fluvial, 24
08015 Parets del Vallès
Barcelona, Spain
Tel: +34 935 710 400



- **Primary labelling**



— The label attached on the device surface


	
Label	Product picture

- **Secondary labelling**


The secondary labelling and its placement in the carton box are shown in the following pictures.

— External label

	
Label	Product picture



— Additional external label

<p>DG Spin</p> <p>6019434</p> <p>3060294</p> <p>GRIFOLS</p> <p>Contenido/Conteúdo: Consulte lista de componentes "packing list"</p> <p>Imported and Marketed by: Grifols India Healthcare Pvt. Ltd. C/o DHL Supply Chain India Pvt. Ltd. A-185, Oldia Industrial Area, Phase - 1, New Delhi, 16, GIGLIANA ROAD, Ftn. 110029 - Customer Care No. 0124-4511620 Email: orders.care@grifols.com Registration No.: Diago-ESP/IND/001291 Product Name: DG Spin - Net Quantity: 1 Unit Size and Dimensions: 525 mm (depth) x 410 mm (width) x 180 mm (height) MOP Rs. 4,75,000/- (inclusive of all taxes)</p> <p>Para uso exclusivo en Laboratorios Clínicos o de Gabinetes. País de origen: España. Importado por: Grifols México, S.A. de C.V. Egecena Caim N. 919 y 913 - Parque Industrial Botones Norte CP. 85150, Zanguay, Jalisco - MEXICO. Distribuido por: Laboratorios LICON, S.A. Avenida Industrial Eléctrica de Mexico N.º 3, Colonia Kocacabaalco, Tlaxiapañita de Bar. CP. 54080 - Estado de Mexico - MEXICO</p> <p>产品名称: 医用离心机 型号: 规格: DG Spin 产品序列号/电源连接条件/生产日期: 参见仪器英文标签 备案号: 国械备20161402号 备案人/生产企业名称: Diagnostic Grifols, S.A. 葛立福诊断股份有限公司 备案人地址/生产企业名称/生产地址: Passeig Fluvial, 24, 08150-Parets del Vallès, Barcelona, Spain 联系方式: (34) 935710400 代理人: 葛立福医药科技(上海)有限公司 联系方式: 021-60201510 代理人住所: 上海市静安区南京西路1515号北楼11层1104、1105、1109- 1104、1105、1106、1107、1108、1109、1110室(实际楼层10层) 使用期限: 5年。但如果按需要对仪器进行正确的维护,并按计划定期 检修,到期后,仪器仍能按其用途安全使用。 其他内容详见说明书</p> <p>Fabricado en España / Made in Spain</p>	
<p>Label</p>	<p>Product picture</p>

KEMENKES RI AKL XXXXXXXXX
 Distributed by:
 PT. Frismed Hoslab Indonesia
 Ruko Bidex Blok G No. 30-31 CBD BSD City,
 Tangerang Selatan - Indonesia 15321
 Phone/Fax : +6221 53167150-52/+6221 5316715

Additional regulatory specific information for Indonesia to be included in the Additional external label (DRAFT)

