

PRODUCT

POWDERED LATEX EXAMINATION GLOVES



DETAILED

TRADEMARK: SUWARNA LATEX EXAMINATION GLOVES

TYPE: • POWDERED

MATERIAL: NATURAL RUBBER LATEX

USAGE:

TO PROTECT PATIENT AND USER FROM CROSS CONTAMINATION
FOR MEDICAL EXAMINATION PURPOSE, THERAPEUTIC, DIAGNOSTIC, TREATMENT, TO HANDLE
CONTAMINATED MATERIAL

ALSO CAN BE USED IN FOOD INDUSTRY, CLEANING, ELECTRONIC INDUSTRY, PHARMACY,
PAINT AND COATING, AGRICULTURE, HOUSEWORK, ETC.

FEATURE:

- ◆ NON STERILE
- ◆ SINGLE USE
- ◆ DISPOSABLE
- ◆ AMBIDEXTROUS
- ◆ MADE OF NATURAL LATEX RUBBER

INNER PACKING: 100 PCS BY WEIGHT / BOX

OUTER PACKING: 10 BOXES / CARTON

QUALITY STANDARD: ASTM D-3578-19

SIZE: X-SMALL, SMALL, MEDIUM, LARGE, X-LARGE

PHYSICAL DIMENSION

DIMENSION	STANDARD	
	SUWARNA	ASTM D3578-19
LENGTH (mm)	230 min.	230 min.
WIDTH (mm)		
X-Small	65 - 75	70 ± 10
Small	75 - 85	80 ± 10
Medium	90 - 100	95 ± 10
Large	105 - 115	111 ± 10
X-Large	115 - 125	120 ± 10
THICKNESS - SINGLE WALL (mm)		
FINGER TIP	0.10 +/- 0.01	0.08 min.
PALM	0.09 +/- 0.01	0.08 min.

PHYSICAL PROPERTIES

PROPERTIES	STANDARD	
	SUWARNA	ASTM D3578-19
ELONGATION AT BREAK (%)	650 min.	650 min.
TENSILE STRENGTH (MPa)	18 min.	18 min.

STORAGE: GLOVES SHALL MAINTAIN THEIR PROPERTIES WHEN STORED IN A DRY PLACE WHICH TEMPERATURE MAY NOT REACH OVER 40°C.
KEEP AWAY FROM DIRECT SUNLIGHT AND OXIDIZING AGENTS.

Shelf-life : 5 years from the date of manufacturing

PRECAUTIONS Disposable Latex Examination Gloves are made for single use only. Safe use of this glove by or on latex sensitized individuals has not been established. Components used in making gloves may cause allergic reactions for some users. If reactions are suspected, cease using the products, wash your hands thoroughly.